

# NanoPump Sanitary High Pressure Pump

Sanitary transfer of small sample of fluid with Pressure

## Unique Benefits

- **Energy efficient:** 1/2 energy loss of the most homogenizers
- **Small:** Integrated design deliver the light weight and small dimensions
- **Smart:** Programming control systems confer diligent functions
- **Silent:** Noiseless performance with high energy efficiency
- **Strong:** More than 30L/hr
- **CE compliant and RoHS compliant**

Piston Pump/Flow controllable



NanoPump series electric Sanitary High Pressure Pump is an equipment operated with the touch screen and controlled intelligently by program. It can be used for laboratory preparation of fat emulsion, liposome, nanosuspension, microemulsion, lipid microsphere, nanoemulsions, dairy products, infusion solutions, cell disruption, juice homogeneity, fine chemical engineering, dye and etc. The maximum working pressure is 30,000 Psi/2000 Bar. All parts touching with medium are FDA approved sanitary materials including 316L stainless steel, 17-4PH stainless steel, titanium alloy, tungsten carbide, PTFE, UHMWPE.

## Specifications

Catalog No.	Nanopump-1500	Nanopump-1000	Nanopump-800	Nanopump-600	Nanopump-300
Flow rate	6L/hr	8L/hr	10L/hr	18L/hr	30L/hr
Min. Sample	20mL	20mL	20mL	50mL	100mL
Max. Pressure	1500bar	1000bar	800bar	600bar	300bar
Dimensions (cm)	85×36×40	85×36×40	85×36×40	85×36×40	85×36×40
Weight	60 kg	60 kg	60 kg	55 kg	50 kg
Max. Temp.	80°C (176°F)				
Power	110V/220V/230V				
Cleaning	Flush to clean/CIP/SIP				
Features	Digital display of pressure or sanitary pressure gauge				



Phone: +1(323) 825-1851  
Fax :+1(323) 978-4693  
Email: Nano@Genizer.com  
[http://www.genizer.com/nanopump\\_p0024.html](http://www.genizer.com/nanopump_p0024.html)

Discount available for public research. Genizer.com reserves the right to change specifications without notice. © 2009 by Genizer. All rights reserved. English Version 091620.

# NanoPump Sanitary High Pressure Pump




NanoPump-1500bar



Control Screen

## 3 Steps for running the NanoGenizer

**Step 2:** Adjust the Speed from Low to High



**Step 3:** Press the Start or Stop

## Standard Features

Parts	1 set of Diamond interaction chamber
Control System	High pressure programming control systems®: Touch Screen, Speed control, Auto stop control by volume, time, pressure or temperature, settable volume control as low as 1mL, display of flow rate and time, overload protection
Pressure Gauge	Electric
Inlet type	1" Triclamp
Outlet type	1" Triclamp
Inlet Reservoir	50mL/100mL/200mL/500mL Syringe; 500mL/1L/2L Stainless Cylinder
Outlet Reservoir	50mL/100mL/200mL/500mL Syringe; 500mL/1L/2L Stainless Cylinder
Product material	316L stainless steel, Tungsten Carbide, Viton, Teflon, UHMWPE
Material standard	Pharmaceutical Grade, FDA & GMP approvable
Warranty	1 year against any manufacturing defects

## Available Option

Parts	High Pressure Extruders, Homogenizing valves, Heat Exchanger
Control System	Mechanical control; Speed, start, stop control
Detector	Pressure gauge, Pressure and Temp. transducer
Outlet type	2" Tri-clamp or Luer
Inlet Reservoir	10ml, 20ml, 50ml, 100ml, 200ml, 500ml Syringe, Tank or Online
Outlet Reservoir	Titanium High pressure Cylinder (Resistant to strong acid and base)



Phone: +1(323) 825-1851  
 Fax : +1(323) 978-4693  
 Email: Nano@Genizer.com  
[http://www.genizer.com/nanopump\\_p0024.html](http://www.genizer.com/nanopump_p0024.html)

Discount available for public research. Genizer.com reserves the right to change specifications without notice. © 2009 by Genizer. All rights reserved. English Version 091620.